





Multi-analyzer Reflectometry

Reflectometry: $k_i = k_f$

$$Q = \frac{4\pi \sin \theta}{\lambda}$$



Non-equilibrium dynamics can be probed on timescales across 3 orders of magnitude.



Operating Principles


- Reflectometry is an elastic scattering technique, $k = k_i = k_f$:
 - Reactor (continuous):** k traditionally determined by monochromating the incident beam.
 - Time of flight (pulsed):** polychromatic beam with k determined by arrival time at the detector.
 - CANDOR:** Combines the advantages of pulsed and continuous sources by using a crystal analyzer to determine k for a polychromatic beam (x20 gain over monochromatic NR).
- Reflectivity curves (R) are collected by rocking the sample to different angles of incidence (θ) and collecting the reflected beam by rotating the detector arm by (2θ).
- Curve fitting the data determines the depth profile of composition, density, and magnetic field with sub-Å precision.

Example: Ultra-low Background


20x speed enhancement over conventional monochromatic reflectometers.

18 orders of magnitude in reflectivity readily achievable, even in difficult sample environments like liquid cells.

Hoogerheide, et. al. JAC 55, 58 (2022)



Example: Polarized Beam



Polarized beam capability on CANDOR enables magnetic measurements with very high precision.

Difference between $\uparrow$ and $\downarrow$ neutron reflectivity clearly distinguishes out to 2 nm in a 12 hr measurement.



CANDOR Neutron Reflectometer



Instrument Webpage: <https://www.nist.gov/ncnr/chrs-candor-white-beam-reflectometer>

CANDOR is the world's fastest neutron reflectometer, with an unprecedented number of neutrons on sample and world-leading signal/noise ratios due to the unique combination of a continuous white beam incident on the sample with wavelength analysis of the scattered beam using cold PG analyzer crystals. CANDOR is ideally suited for studying the time-resolved nuclear and magnetic structures of thin films important for fields such as advanced computing, quantum communications and encryption, storage technologies, and biological functions across membranes.

Beam Preparation

Chromatic Analysis

Be/Bi Filter

8" Be: eliminates < 4 Å neutrons
8" Bi: eliminates gammas

Shutter

Guide
60mm x 150mm

4.75 Å
Monochromator

Converging guide
increasing neutrons on
sample using tunable
beam divergence


RF
Flipper

Broad energy scattered
(from the sample)
neutron beam

> 6 Å Deflector


High-efficiency
double-V polarizer:
 $P > 99\%$

Sample
Stage



Scattered beam is
diffracted into thin film
 $^6\text{LiF}\cdot\text{ZnS}$ (Ag) neutron
detectors invented
specifically for CANDOR

Detector banks are
cooled below 20 K to
minimize background



Cryostats

System

- Mini-cryostat
- Furnace
- Cryostat-compatible multi-sample holder

Range

4 K - 320 K
200 °C - 600 °C
4 positions

Pressure and Humidity

System

- RH and temp. control
- Gas Loading
- Pressure Cells

Range

up to 100% and 90 °C
H, D, CH₄, CO₂
1 GPa He, Ne, Ar

Magnets

System

- Permanent magnet mount
- Electromagnet
- Superconducting magnet
- Helmholtz coils

Range

20 mT
700 mT
3T

Other

- Sample changing and liquid handling for solid/liquid interfaces.
- In-situ voltage application up to 5 kV

Operational Range

- **White Beam:** wavelengths from 4 Å - 6 Å
- **Monochromatic Beam:** 4-bounce monochromator fixed at 4.75 Å
- **Q range:** 0.05 nm⁻¹ - 30 nm⁻¹

Beam Divergence

- Two variable incident slits, with background reduction slit for scattered beam.
- Converging guide with variable beam divergence.
- Two operational detector banks, with option for 18 more covering 2° of scattered beam each.

Polarization

- Before sample: Double-V supermirror polarizer and polarized guide path.
- After sample: ³He polarization analysis
- RF flippers and guide fields