

SOP No. 6

Recommended Standard Operations Procedure  
for  
Weighing by Transposition

1. Introduction

1.1. Purpose

This method eliminates possible error in comparison of weights when using an equal-arm balance, due to a small inequality in the length of the arms. While this procedure is written to describe the measurement of masses of weights for calibration purposes, it may be used for general mass determinations, whenever accurate weighing is required. (Note: This is a single - transposition procedure, hence does not eliminate the effects of linear drift).

1.2. Prerequisites

- 1.2.1. Verify that standard weights are available with valid certificates of calibration.
- 1.2.2. Verify that the balance to be used is in good working order. This may be verified by maintenance of a valid control chart. Otherwise, a series of trial measurements should be made to ascertain that the balance is functional and that its sensitivity is adequate for the requirements of the measurement.
- 1.2.3. Verify the ability of the technician to make precision weighings and that he/she has had training in this specific procedure.

2. Methodology

2.1. Scope, Precision, Accuracy

The measurement principle is applicable to weighing with any equal-arm balance. The range of weights that can be calibrated will depend upon the capacity of the balance used. The precision of intercomparison will depend upon the sensitivity of the balance, and the care exercised in making the required weighings. The accuracy will depend on the accuracy of calibration of the standard weights and the precision of intercomparison.

2.2. Summary

The masses of a standard and test weight are intercompared. The loads are then interchanged and a second intercomparison is made. A small calibrated weight is then added to one of the masses and a

third weighing is made to determine the sensitivity of the balance under the load conditions.

### 2.3. Apparatus/Equipment

- 2.3.1. Equal arm balance of capacity and sensitivity sufficient for the weights tested. The index scale of the balance is conveniently numbered 0 to 20 with 10 as the center division, although other numbering systems such as 0 to 200 are possible. The graduations should be numbered so that addition of a weight to the left arm will increase the balance readings. A system in which the center division is 0 is not recommended since the negative readings that result can cause observational and computational problems.
- 2.3.2. Standard weights with valid calibrations, traceable to NBS.
- 2.3.3. Small calibrated weights (usually decimal fractions) to be used as tare weights.
- 2.3.4. Equipment capable of loading and unloading weights on balance without damage to either (especially important in the case of large weights).

### 2.4. Symbols

The following symbols are used in this procedure.

- S - standard weight
- X - weight calibrated
- T - counterweight
- t - small calibrated. A subscript s or x is used to indicate the larger weight with which it is associated.
- sw - small calibrated weight used to evaluate the sensitivity of the balance.
- M - the mass of a specific weight. Subscripts s, x, t, sw, are used to identify the weight.
- AM - the apparent mass of a specific weight. Subscripts s, x, t, sw together with numerical sub-subscripts as necessary are used to identify the weight.

## 2.5. Procedure

2.5.1. The following weighing sequences may be used:

### Option A

Measurement No.	Left Arm	Right Arm	Observation
1	$X + t_{x_1}$	$S + t_{s_1}$	$O_1$
2	$S + t_{s_2}$	$X + t_{x_2}$	$O_2$
3	$S + t_{s_2} + sw$	$X + t_{x_2}$	$O_3$

### Option B

Measurement No.	Left Arm	Right Arm	Observation
1	$S + t_{s_1}$	$X + t_{x_1}$	$O_1$
2	$X + t_{x_2}$	$S + t_{s_2}$	$O_2$
3	$X + t_{x_2} + sw$	$S + t_{s_2}$	$O_3$

## 2.5.2. Preliminary Procedures

Conduct preliminary measurements (without recording data) to determine the values for  $t_s$  and  $t_x$  as described in 2.5.3. This will serve to warm up the balance and to facilitate the measurements. In a series of calibrations, this step may be minimized or eliminated after the first series of measurements. Depending on the corrections for S and X, weights  $t_s$  and  $t_x$  may be unnecessary in some cases.

## 2.5.3. Measurement Procedure (Option A)

All observations should be recorded on suitable data sheets such as those in the Appendix. Record the laboratory ambient temperature, barometric pressure, and relative humidity.

2.5.3.1. Observation 1. Place the test weight(s) on the left arm and the standard(s) on the right arm. Add  $t_x$ , and/or  $t_s$  (small calibrated weights) as necessary to obtain an approximate balance. Read and record all data and the sum of the turning points  $O_1$  (see GMP No. 1).

2.5.3.2. Observation 2. Place the standard weight(s) on the left arm and the test weight(s) on the right arm with appropriate small weights, as was done in Observation 1. The tare weights carried with S and X may have to be changed so that  $O_2$  is within one division of  $O_1$  on a 0 to 20 division graduated

scale. Read and record all data and the sum of the turning points,  $O_2$ .

2.5.3.3. Observation 3. Without disturbing weights already in place from Observation 2, add a small calibrated sensitivity weight ( $sw$ ) to the left arm, sufficient to displace the turning points 3 to 5 divisions on a 20 division graduated scale. Record the sum of the turning points,  $O_3$ .

2.5.4. Measurement Procedure (Option B)

Conduct 3 measurements as in 2.5.3.1, 2.5.3.2, and 2.5.3.3 using the weighing schedule Option B.

3. Calculations

3.1. No Air Buoyancy Correction

The following equations may be used to calculate the corrections required for the test weight(s). The calibrated apparent mass values for  $t_s$ ,  $t_x$ , and  $sw$  must be used.  $C_s$  is the apparent mass correction for the standard weight(s),  $S$ . The symbols  $N_s$  and  $N_x$  refer to the nominal values of  $S$  and  $X$ , respectively.

Option A

$$C_x = C_s + \frac{AM_{t_{s1}} + AM_{t_{s2}} - AM_{t_{x1}} - AM_{t_{x2}}}{2} + \frac{O_1 - O_2}{2} \cdot \frac{AM_{sw}}{O_3 - O_2} + N_s - N_x$$

Option B

$$C_x = C_s + \frac{AM_{t_{s1}} + AM_{t_{s2}} - AM_{t_{x1}} - AM_{t_{x2}}}{2} + \frac{O_2 - O_1}{2} \cdot \frac{AM_{sw}}{O_3 - O_2} + N_s - N_x$$

3.2. Air Buoyancy Correction

3.2.1. Calculate the air density,  $\rho_A$ , as described in Section 8 of the Appendix to SOP No. 2 or obtain from Table 9.9.

3.2.2. Calculate the mass,  $M_x$ , of the test weight  $X$ , using the mass of the standard weight(s), the tare weights, and the sensitivity weights according to the optional sequence used. It is assumed that the densities of the combination of tare weights used do not differ significantly in the weighing design.

Symbols used:  $\rho_A$  = air density  
 $M_i$  = mass for weight  $i$   
 $\rho_i$  = density for weight  $i$

3.2.2.1. Optional Sequence A

$$M_x = \frac{M_s \left(1 - \frac{\rho_A}{\rho_s}\right) + \frac{O_1 - O_2}{2} \cdot \frac{M_{sw} \left(1 - \frac{\rho_A}{\rho_{sw}}\right)}{O_3 - O_2} + \frac{M_{ts1} + M_{ts2}}{2} \left(1 - \frac{\rho_A}{\rho_{ts}}\right) - \frac{M_{tx1} + M_{tx2}}{2} \left(1 - \frac{\rho_A}{\rho_{tx}}\right)}{1 - \frac{\rho_A}{\rho_x}}$$

3.2.2.2. Optional Sequence B

$$M_x = \frac{M_s \left(1 - \frac{\rho_A}{\rho_s}\right) + \frac{O_2 - O_1}{2} \cdot \frac{M_{sw} \left(1 - \frac{\rho_A}{\rho_{sw}}\right)}{O_3 - O_2} + \frac{M_{ts1} + M_{ts2}}{2} \left(1 - \frac{\rho_A}{\rho_{ts}}\right) - \frac{M_{tx1} + M_{tx2}}{2} \left(1 - \frac{\rho_A}{\rho_{tx}}\right)}{1 - \frac{\rho_A}{\rho_x}}$$

3.2.2.3. Calculate the mass correction,  $C'_x$ , as follows:

$$C'_x = M_x - N_x$$

when  $N_x$  is the nominal value for X

3.2.2.4. Calculate the apparent mass of X versus the desired reference density of 8.0 g/cm<sup>3</sup> or brass. It is recommended that the apparent mass versus 8.0 g/cm<sup>3</sup> be reported unless otherwise requested. The density of X,  $\rho_x$ , must be in g/cm<sup>3</sup>.

3.2.2.4.1. Apparent mass versus 8.0 g/cm<sup>3</sup>

$$AM_x \text{ vs. } 8.0 = M_x \frac{\left(1 - \frac{.0012}{\rho_x}\right)}{0.999850}$$

3.2.2.4.2. Apparent mass versus brass

$$AM_x \text{ vs. brass} = M_x \frac{\left(1 - \frac{.0012}{\rho_x}\right)}{0.999857}$$

#### 4. Assignment of Uncertainty

The limits of uncertainty,  $U$ , include estimates of the uncertainty of the mass standards,  $U_s$ , plus the uncertainty of measurement,  $U_m$ , at the 99.73% level of confidence. The latter is estimated by

$$ts$$

where  $s$  is the standard deviation of measurement and  $t$  is obtained from Table 9.3.

Then

$$U = \pm [U_s + ts]$$

##### 4.1. Precision of Measurement Known from Control Chart Performance (SOP No. 9)

The value for  $s$  is obtained from the control chart data for transposition measurements. Statistical control will need to be verified by measurement of at least one check standard while the above measurements are in progress.

Use the value of  $t$  (corresponding to a probability level of 99.73%) from Table 9.3 appropriate for the number of degrees of freedom,  $\nu$ , on which the control limits of the control chart are based.

##### 4.2. Precision Estimated from a Series of Measurements

Measure a stable test object at least 7 times, no two measurements of which may be made on a single day. Calculate the mean and a standard deviation in the conventional manner. The latter is the value of  $s$  that is to be substituted in the equation given in Section 4. In this case select the value for  $t$  from Table 9.3 based on the number of degrees of freedom involved in computing  $s$ .

Note: Repetitive measurements made on the same day evaluate only the short-term standard deviation.

#### 5. Report

##### 5.1. Report results as described in SOP No. 1, Preparation of Calibration/Test Reports

Appendix

**Transposition Data Sheet  
for Equal-Arm Balance (Option A)**

Test No.: \_\_\_\_\_ Sheet No.: \_\_\_\_\_ Date: \_\_\_\_\_

Item Identification: \_\_\_\_\_ Balance: \_\_\_\_\_

Standard Identification: \_\_\_\_\_ Observer: \_\_\_\_\_

Temperature: \_\_\_\_\_ Pressure: \_\_\_\_\_ Rel. Hum.: \_\_\_\_\_

$C_s(C'_s)$ - \_\_\_\_\_  $\pm$  \_\_\_\_\_  $AM_{sw}(M_{sw})$ - \_\_\_\_\_

$\rho_s$ - \_\_\_\_\_  $\rho_x$ - \_\_\_\_\_  $\rho_{sw}$ - \_\_\_\_\_

Time \_\_\_\_\_ Balance standard deviation = \_\_\_\_\_

Measure- ment No.	Weights		Turning Points		Sum
	Left Arm	Right Arm	Low	High	
1	$X + t_{x1}$ $t_{x1}$ -	$S + t_{s1}$ $t_{s1}$ -			$O_1$ -
2	$S + t_{s2}$ $t_{s2}$ -	$X + t_{x2}$ $t_{x2}$ -			$O_2$ -
3	$S + t_{s2} + sw$	$X + t_{x2}$			$O_3$ -

Time \_\_\_\_\_

$AM_{t_{s1}}(M_{t_{s1}})$ - \_\_\_\_\_  $\pm$  \_\_\_\_\_

$AM_{t_{x1}}(M_{t_{x1}})$ - \_\_\_\_\_  $\pm$  \_\_\_\_\_

$AM_{t_{s2}}(M_{t_{s2}})$ - \_\_\_\_\_  $\pm$  \_\_\_\_\_

$AM_{t_{x2}}(M_{t_{x2}})$ - \_\_\_\_\_  $\pm$  \_\_\_\_\_

$\rho_{t_s}$ - \_\_\_\_\_

$\rho_{t_x}$ - \_\_\_\_\_

Appendix

Transposition Data Sheet  
for Equal-Arm Balance (Option B)

Test No.: \_\_\_\_\_ Sheet No.: \_\_\_\_\_ Date: \_\_\_\_\_

Item Identification: \_\_\_\_\_ Balance: \_\_\_\_\_

Standard Identification: \_\_\_\_\_ Observer: \_\_\_\_\_

Temperature: \_\_\_\_\_ Pressure: \_\_\_\_\_ Rel. Hum.: \_\_\_\_\_

$C_s(C'_s) =$  \_\_\_\_\_  $\pm$  \_\_\_\_\_  $AM_{sw}(M_{sw}) =$  \_\_\_\_\_

$\rho_s =$  \_\_\_\_\_  $\rho_x =$  \_\_\_\_\_  $\rho_{sw} =$  \_\_\_\_\_

Time \_\_\_\_\_ Balance standard deviation = \_\_\_\_\_

Measure- ment No.	Weights		Turning Points		Sum
	Left Arm	Right Arm	Low	High	
1	$S + t_{s1}$ $t_{s1}^-$	$X + t_{x1}$ $t_{x1}^-$			$O_1^-$
2	$X + t_{x2}$ $t_{x2}^-$	$S + t_{s2}$ $t_{s2}^-$			$O_2^-$
3	$X + t_{x2} + sw$	$S + t_{s2}$			$O_3^-$

Time \_\_\_\_\_

$AM_{t_{s1}}(M_{t_{s1}}) =$  \_\_\_\_\_  $\pm$  \_\_\_\_\_

$AM_{t_{x1}}(M_{t_{x1}}) =$  \_\_\_\_\_  $\pm$  \_\_\_\_\_

$AM_{t_{s2}}(M_{t_{s2}}) =$  \_\_\_\_\_  $\pm$  \_\_\_\_\_

$AM_{t_{x2}}(M_{t_{x2}}) =$  \_\_\_\_\_  $\pm$  \_\_\_\_\_

$\rho_{t_s} =$  \_\_\_\_\_

$\rho_{t_x} =$  \_\_\_\_\_



Appendix

Transposition Data Sheet  
for Equal-Arm Balance (Option A)

Test No.: 135 Sheet No.: 1 Date: 8/26/86

Item Identification: 1000 lb No. 4211 Balance: Russell

Standard Identification: #11 and #22 Observer: HO

Temperature: 23.5°C Pressure: 745.8 mm Hg Rel. Hum.: 50%

$C_s(C'_s) =$ 0.0047 lb  $\pm$  0.0030 lb  $AM_{sw}(M_{sw}) =$ 0.01 lb

$\rho_s =$ 8.0 g/cm<sup>3</sup>  $\rho_x =$ 7.0 g/cm<sup>3</sup>  $\rho_{sw} =$ 7.8 g/cm<sup>3</sup>

Time 1:15 p Balance standard deviation = 0.0018 lb

Measurement No.	Weights		Turning Points		Sum
	Left Arm	Right Arm	Low	High	
1	$X + t_{x1}$	$S + t_{s1}$	4.0	<del>16.3</del>	$O_1 = 20.35$
	$t_{x1} = 3.6171\text{lb}$	$t_{s1} = 3.6251\text{lb}$	4.1	16.3	
2	$S + t_{s2}$	$X + t_{x2}$	<del>5.6</del>	<del>14.1</del>	$O_2 = 19.75$
	$t_{s2} =$	$t_{x2} = 3.6121\text{lb}$	5.7	14.1	
3	$S + t_{s2} + sw$	$X + t_{x2}$	<del>8.4</del>	<del>16.0</del>	$O_3 = 24.3$
			8.4	15.9	
				15.9	

Time 1:40 p

$AM_{t_{s1}} (M_{t_{s1}}) =$ 3.625  $\pm$  .0000751b  $AM_{t_{x1}} (M_{t_{x1}}) =$ 3.617  $\pm$  negligible

$AM_{t_{s2}} (M_{t_{s2}}) =$ 3.625  $\pm$  .0000751b  $AM_{t_{x2}} (M_{t_{x2}}) =$ 3.612  $\pm$  negligible

$\rho_{t_s} =$ 8.0 g/cm<sup>3</sup>

$\rho_{t_x} =$ 7.8 g/cm<sup>3</sup>

$$C_x = C_s + \frac{AM_{t_{s1}} + AM_{t_{s2}} - AM_{t_{x1}} - AM_{t_{x2}}}{2} + \frac{O_1 - O_2}{2} \cdot \left( \frac{AM_{sw}}{O_3 - O_2} \right)$$

$$C_x = .0047 + \frac{3.625 + 3.625 - 3.617 - 3.612}{2} + \frac{20.35 - 19.75}{2} \cdot \frac{0.01}{24.3 - 19.75}$$

$$C_x = .0047 + .0105 + .000659$$

$$C_x = 0.015859 \text{ lb}$$

Uncertainty =  $t_s + U_s$                       Degrees of freedom in  $s = 20$

$$\text{Uncertainty} = 3.422 (0.0018) + (.0030 + .000075)$$

$$\text{Uncertainty} = 0.0092346 \text{ lb}$$

$$C_x = 0.0159 \pm 0.0092 \text{ lb}$$