

NIST DNA Analyst Webinar Series


Validation and use of likeLTD

Matthew Greenhalgh
Director of Forensic Science

18 September 2014


Twenty Five Years of DNA Analysis

- Established ICI (Astrazeneca) in 1987
- Specialist Forensic DNA analysis
- One of Europe's largest paternity testers
- Contracted to >80% of UK police forces
- Approx 475 UK employees
- A LabCorp company since 2011



Background to the problem


- The increasing sensitivity of DNA analysis methods and their use in a wider range of case types has resulted in more mixtures:



Interpretation procedures prior to likeLTD


- Match probability for single source profiles (or major/minor)
- LR method for 2 person mixtures where both individuals are fully represented
- “Dlugosz” expert opinions without statistical evaluation of the match in a limited number of cases
 - Court of Appeal Ruling R v Dlugosz, R v Pickering and R v MDS ([2013] EWCA Crim 2)

Distribution of DNA results - volume crime




Based on 9165 samples (3 months data)

Breakdown of mixture results




Profiles suitable for analysis



- Low level mixture
- 2 contributors?
- Suspect + unknown
- Allele dropout
- 3 replicates

Profiles suitable for analysis (*continued*)



- 3 contributors?
- Suspect + victim + unknown
- No apparent dropout
- 2 replicates

DNA Resolve – Cellmark Forensic Services

- Uses likeLTD software written by Prof. Balding (UCL London) combined with a Cellmark user-interface
- Models allele drop-out using Tvedebrink statistical model
- Option to include allele drop-in
- Can deal with a maximum of 2 unprofiled contributors in a mixture
- Multiple replicates can be analysed

likeLTD input files (.csv)

- Crime stain file


| Stain | Profiling system | Plate/Run | | D3 | vWA | D16 | D2 | D8 |
|-------|------------------|-----------|-----------|----------------|-------------------|-------------|-------------|----------------|
| CSP | SGM+ | F2 | Allelic | 14,15,16,17,18 | 14,15,16,17,18 | 10,11,12,13 | 17,18,24 | 10,13,14,15 |
| CSP | SGM+ | | Uncertain | 13 | | | 16 | |
| CSP | SGM+ | G2 | Allelic | 14,16 | 15,16,17,18 | 10,11,12,13 | 17,18,25 | 10,12,13,14,15 |
| CSP | SGM+ | | Uncertain | 15 | | 9 | 16,19 | |
| CSP | SGM+ | H2 | Allelic | 14,15,16,17 | 14,16 | 10,11,12,13 | 17,18 | 10,13,14,15 |
| CSP | SGM+ | | Uncertain | | 15 | | 25 | |
| CSP | SGM+ | B3 | Allelic | 14,15,16,17,18 | 14,15,16,17,18,19 | 9,10,11,13 | 17,18,20,25 | 10,12,13,14,15 |
| CSP | SGM+ | | Uncertain | 13 | | 12 | 24 | |

- Reference File

| Individual | known/queried | D3 | vWA | D16 | D2 | D8 | D21 |
|------------|---------------|-------|-------|-------|-------|-------|---------|
| Suspect | queried | 14,16 | 14,17 | 11,13 | 18,24 | 12,14 | 28,31.2 |
| Victim | known | 15 | 16,20 | 9,11 | 23 | 12,13 | 28,33.2 |


Allelic and Uncertain peaks

- Allelic Peaks



- Uncertain Peaks

- Sub-threshold
- Possible Stutter



Allele report

Data provided by forensic scientist

Crime scene profiles (CSP)


| ru n | status | D3 | vWA | D16 | D2 | D8 | D21 | D18 | D19 | TH01 | FGA |
|---------|-----------|----------------|-------------------|-------------|-------------|----------------|--------------------|----------------|---------------|-----------|----------------|
| 1 | certain | 14,15,16,17,18 | 14,15,16,17,18 | 10,11,12,13 | 17,18,24 | 10,13,14,15 | 27,28,32.2 | 12,13,14,17 | 13,14 | 6,7,9,9.3 | 21,23,25 |
| - | uncertain | 13 | | | 16 | | | 11,16 | 12 | | 20,26 |
| 2 | certain | 14,16 | 15,16,17,18 | 10,11,12,13 | 17,18,25 | 10,12,13,14,15 | 28,30,31,31.2,32.2 | 16,17 | 13,14,15 | 6,7,9 | 19,21,22,23 |
| - | uncertain | 15 | | 9 | 16,19 | | 27 | | | | |
| 3 | certain | 14,15,16,17 | 14,16 | 10,11,12,13 | 17,18 | 10,13,14,15 | 28,32.2 | 12,14,15,17,18 | 13,14,14.2 | 6,7,9,9.3 | 21,22,23 |
| - | uncertain | | 15 | | 25 | | | 16 | | | 20,25,26 |
| 4 | certain | 14,15,16,17,18 | 14,15,16,17,18,19 | 9,10,11,13 | 17,18,20,25 | 10,12,13,14,15 | 27,28,32.2 | 13,16,17 | 13,14,15,15.2 | 6,7,9,9.3 | 20,21,22,23,25 |
| - | uncertain | 13 | | 12 | 24 | | | | | | 24 |

Reference profiles

| profile | D3 | vWA | D16 | D2 | D8 | D21 | D18 | D19 | TH01 | FGA |
|---------|-------|-------|-------|-------|-------|---------|-------------|--------------|------|-------|
| Q | 14,16 | 15,16 | 10,13 | 17,18 | 10,15 | 28,32.2 | 17 | 14 | 6,7 | 21,23 |
| K | 15 | 14,18 | 11,9 | 25,24 | 14 | 30,31.2 | 13,15 | 14,15 | 7,9 | 22,25 |
| Other | 17,18 | 17,19 | 12 | 20 | 13,12 | 27,31 | 12,14,16,18 | 13,14.2,15.2 | 9.3 | 19,20 |

Alleles that are replicated, unreplicated or absent in the crime scene profile, using the certain designations only.

Allele report (*continued*)



The number of 'certain' alleles that cannot be attributed to the known profile(s).

Unusual alleles

| source | locus | allele | EA1.freq | EA3.freq | EA4.freq | error |
|-----------------------|-------|--------|----------|----------|----------|-------|
| Crime scene uncertain | D3 | 13 | 5 | 7 | 0 | - |
| Crime scene uncertain | D18 | 11 | 10 | 1 | 6 | - |

Alleles are automatically checked against the database. An error will be reported if an allele is absent from the database, or present more than once, or if a locus is absent.

Approximate representation

| Contributor | Rep 1 | Rep 2 | Rep 3 | Rep 4 | Total |
|-------------|-------|-------|-------|-------|-------|
| Q | 100 | 100 | 94 | 100 | 99 |
| K | 61 | 61 | 50 | 78 | 62 |

Allele report (*continued*)

The fraction of an individual's alleles (as a percentage) that have been designated as 'certain' alleles in each replicate. This estimate is not used by likeLTD, and is intended to assist informal assessments of possible known contributors to the CSP. A more formal approach is to do a likeLTD run to compute the likelihood ratio (LR) for that individual contributor.

Suggested parameter values

| nUnknowns | doDropin | Recommendation |
|-----------|----------|----------------|
| 2 | FALSE | recommended |

Recommended values for 'nUnknowns', choose from 0, 1 or 2 (likeLTD automatically adds an additional unknown X to the defence hypothesis in place of the queried profile Q).

Recommended values for 'doDropin', choose from 'TRUE' or 'FALSE'.

All the attributable alleles must either come from an unknown or dropin.

System information

| Type | Details |
|------------------------|--|
| Date report generated: | Fri Aug 01 09:07:22 2014 |
| Package | likeLTD |
| Title | Tools to determine DNA profile evidence. |
| Description | Tools to determine DNA profile Weight of Evidence. For further information see the likeLTD guide at the URL provided, or the paper under citation. |
| Depends | R (≥ 2.10), DEoptim, ggplot2, gtools, rtf |
| Suggests | svUnit, scales |
| Imports | gdata, tools, tcltk |
| Version | 5.4.0 |
| Date | 2013-03-15 |
| Author | David Balding, Adrian Timpson, Christopher Steele, Mayeul d'Avezac, James Hetherington. |
| Maintainer | Christopher Steele <c.steele.11@ucl.ac.uk> |
| License | GPL-3 |
| URL | https://sites.google.com/site/baldingstatisticalgenetics/ |

Evaluation report

CFS-920055-14-ver5.4-Evaluation-Report

CFS-920055-14-ver5.4

Prosecution Hypothesis: Q (Q) + K + 2U

Defence Hypothesis: Unknown (X) + K + 2U

Evaluation report (*continued*)

The fraction of an individual's alleles (as a percentage) that have been designated as 'certain' alleles in each replicate. This estimate is not used by likeLTD, and is intended to assist informal assessments of possible known contributors to the CSP. A more formal approach is to do a likeLTD run to compute the likelihood ratio (LR) for that individual contributor.

Likelihoods at each locus

| Likelihood | D3 | vWA | D16 | D2 | D8 | D21 | D18 | D19 | TH01 | FGA |
|-------------------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
| Prosecution.log10 | -2.443 | -4.427 | -1.216 | -3.158 | -1.751 | -6.926 | -9.377 | -6.054 | -1.503 | -4.713 |
| Defence.log10 | -3.432 | -5.615 | -2.404 | -4.444 | -2.993 | -8.446 | -9.455 | -5.790 | -2.039 | -5.821 |
| Ratio.log10 | 0.989 | 1.187 | 1.188 | 1.286 | 1.241 | 1.520 | 0.078 | -0.265 | 0.537 | 1.108 |
| Ratio | 9.739 | 15.397 | 15.407 | 19.324 | 17.432 | 33.120 | 1.196 | 0.543 | 3.441 | 12.827 |

Overall Likelihood

| calculation | estimate |
|-------------------|-----------|
| Prosecution.log10 | -41.568 |
| Defence.log10 | -50.437 |
| Ratio.log10 | 8.869 |
| Ratio | 739308745 |

Theoretical maximum LR

| calculation | estimate |
|------------------------|----------------|
| likelihood ratio | 10915470994030 |
| Log10 likelihood ratio | 13.038 |

Dropout and degradation parameter estimates

Input required using R interface

```
> require(likeLTD)
Loading required package: likeLTD
Loading required package: DEoptim

DEoptim package
Differential Evolution algorithm in R
Authors: D. Ardia, K. Mullen, B. Peterson and J. Ulrich

Loading required package: ggplot2
Loading required package: gtools
Loading required package: rtf
> setwd("P:/Projects/Active/Complex DNA profile interpretation using likeLTD software/")
> datapath <- "."
> admin <- pack.admin.input(
+   cspFile = "CSP1.csv",
+   refFile = "References.csv",
+   databaseFile = "Cellmark_NDNAD-allele-freqs-DNA17_v2_(SGM PLUS ONLY)-wbpl.txt",
+   caseName = "CSP1 batch test"
+ )
> allele.report(admin)
>
> args <- list(
+   nUnknowns = 0,
+   doDropin = FALSE,
+   ethnic = "EA1",
+   adj = 1,
+   fst = 0.02,
+   relatedness = c(0,0)
+ )
>
> hypP <- do.call(prosecution.hypothesis, append(admin,args))
> hypD <- do.call(defence.hypothesis, append(admin,args))
>
> paramsP <- optimisation.params(hypP)
> paramsD <- optimisation.params(hypD)
>
> results <- evaluate(paramsP, paramsD, progBar = FALSE, interim = FALSE)
```

Cellmark user interface

7 P:/Projects/Active/Complex DNA profile interpretation using like.LTD ...

File

Cellmark
Forensic Services

Copyright (C) 2013, Orchid Cellmark,
Abingdon Business Park, Abingdon, Oxon OX14 1DY
This software is developed as a PRIVATE COPY
for Cellmark under the terms of the GNU GPL
NOT FOR DISTRIBUTION OUTSIDE CELLMARK

Set Directory

Case name: CSP7

Frequency database: P:/Projects/Active/Complex ... EA1

File of CSP alleles: CSP7.csv

File of Reference profiles: References.csv

Generate initial allele report

Unprofiled contribs (PH)

Number of unprofiled contributors 1

IBD (1) IBD (2)

Relationship (DH)
(5 to first known contributor) 0 0

Dropins: ☐

Coancestry coefficient (Fst): 0.02

Seed: 42

Calculate Exit

STATUS: Ready

likeLTD validation

- Establishing that GUI did not affect results
- Repeat tests published by Balding
- Additional testing
- Validated under ISO 9001 certification
- Planning to add to ISO 17025 scope

Tests from likeLTD guide

- CSP1** Full profile match to reference Q, single contributor
- CSP2** The two replicates of CSP2 differ from reference Q due to 1 drop in and 2 dropouts
- CSP3** One further drop in and two more dropouts have been introduced
- CSP4** Two contributors: All the alleles of both contributors present in both replicates with no drop in or dropout
- CSP5** Introduces random 50% dropout for the alleles of unknown 1 (U1) not shared with Q

Tests from likeLTD guide (*continued*)

- CSP6** The opposite situation is considered where there is 50% dropout of the alleles of Q not shared with U1
- CSP7** In addition to the 50% dropout for the alleles of Q, 50% of the alleles of U1 generate stutter peaks that are classified as uncertain
- CSP8** Random 50% dropout affects both the alleles of Q and U1


Additional validation tests

- Effects on the LR value of increasing levels of F_{st} for Caucasian, Afro Caribbean and Asian frequency databases + comparison against the reciprocal of the match probability calculated using Cellmark's in house software
- Effect of analysis assuming the defence scenario that the donor of the DNA is the suspect's brother
- Effect of using an incorrect number of contributors

Additional validation tests (*continued*)


- Comparison of likeLTD against Mixture Analyser software
- Comparison of likeLTD against another probabilistic software package (STRmix)
- Effects of using the “uncertain” option for allele calls and varying the number of replicates
- Random reference profiles compared against crime scene stain
 - Single contributor
 - Two contributors

Effect of increasing the number of uncertain and negative loci




Effect of increasing numbers of replicates


**Reference Q vs. CSP4 with 3 neg/uncertain alleles calls.
Effect of increasing numbers of replicates**



Two person mixtures with reference profile included



Two person mixtures with reference profile excluded



Advantages of likeLTD

- Provides an objective LR value in complex cases
- Can incorporate replicate PCR runs
- Freely available open source software
- Theoretical aspects published in peer reviewed journals
- Evidence successfully presented in court

Casework issues

- Limited to a maximum of 2 unknown contributors under Hp
- Memory requirements increase significantly with more loci
- Complex profiles can take several days to analyse on a standard desktop PC
- Can be difficult for non “R” code users

likeLTD publications




Interpreting low template DNA profiles

David J. Balding^a, John Buckleton^b

^a Department of Epidemiology and Public Health, Imperial College, St Mary's Campus, Norfolk Place, London W2 1PG, UK

^b ESR Private Bag 92021, Auckland, New Zealand




Verifying likelihoods for low template DNA profiles using multiple replicates

Christopher D. Steele^{a,*}, Matthew Greenhalgh^b, David J. Balding^a

^a UCL Genetics Institute, Darwin Building, Gower Street, London WC1E 6BT, UK

^b Orchid Cellmark Ltd., Abingdon Business Park, Blacklands Way, Abingdon OX14 1YX, UK



The likeLTD software: an illustrative analysis, explanation of the model, results of performance tests and version history

David J. Balding, Chris Steele
UCL Genetics Institute
Darwin Building, Gower Street
London WC1E 6BT
d.balding@ucl.ac.uk

March 26, 2014

This version coincides with Release 5-3, which improves the allele and output reports, which now produce .doc files, with the same functionality as previous reports. These include improved layout and presentation, and modifications to the allele report. See Sections 1.2 & 1.4.

The previous version (Release 5-2) introduced the function `get.likely.genotypes` that returns the most probable genotypes for each locus, and the most probable whole-profile genotype. See Section 1.5 for more information.

Abstract

likeLTD in court

- likeLTD evidence has been accepted without challenge in more than 10 trials in the UK
- There have been several admissibility challenges (voir dire) - all rejected
 - Evidence was originally ruled inadmissible in the case of R v MDS. At the subsequent retrial the evidence was again challenged but accepted

Summary

- likeLTD has been introduced into forensic casework following internal validation
- Provides objective LR values in complex mixture cases
- Evidence has been accepted in UK courts

Currently evaluating fully continuous probabilistic software

NIST DNA Analyst Webinar Series

Validation and use of likeLTD

Matthew Greenhalgh
Director of Forensic Science

18 September 2014