




Color MRI contrast agents using magnetic microstructures

G. Zabow

S. Dodd, J. Moreland, A. Koretsky

NIH & NIST



U.S. Department of Health and Human Services
National Institutes of Health






Magnetic Resonance Imaging (MRI)









- MRI equipment market ~ \$ 5 Billion / yr (and growing...)
- 30 million MRI exams in the U.S. in 2008 alone.

MRI contrast agents





Topics in Magnetic Resonance Imaging 7(3):181-195 © 1995 Raven Press, Ltd., New York


Contrast agents >> \$ 1 Billion / yr






From Grayscale Contrast to Color

- Color gives extra dimension
- Different regions "tagged" with different colors
- Multiplexing / increased throughput
- Imagine: cancerous / diseased tissue appearing red while normal tissue appears blue....









- Massively parallel wafer-level microfabrication
- Shape determines color.
- Can be turned on or off, sensing physiological conditions
- Structure composition can be from same (FDA approved) materials used in regular contrast agents

U.S. Department of Health and Human Services National Institutes of Health





Commercial Applications

- Provides unique color contrast, and can also double as regular MRI contrast agent (similar to, e.g. Feridex)
- On/off functionality allows for in-vivo real-time physiological sensing, i.e. dynamic contrast.
- Agent sensitivity can exceed that of regular gadolinium MRI agents
- Fluid flow monitoring / stent applications
- Microscale RFID tags, RFID-based microfluidics





Technology stage

- Early stage technology
- Patent pending
- But, already substantial interest:
 - Initial Nature paper, broad media coverage:







- Selected by Nature editors as one of their favorite papers of 2008.
- Follow-up paper selected as cover article in *Nanotechnology* (Sept. 2009).





Collaboration Opportunities

Do you have expertise in any of these fields?

- MRI contrast agents
- MRI calibration
- Micro-RFID tags and microfluidics
- Stent technology
- Flow Cytometry
- MEMS Technology
- Detection/labeling/tracking of biologics

Then, let's talk....





Contact Information

For further information contact:

Licensing:

John Stansberry NIH Office of Technology Transfer stansbej@mail.nih.gov 301-435-5236

Collaboration:

Melissa Maderia, Ph.D.
NCI Technology Transfer Center, CSC for NINDS
maderiam@mail.nih.gov
301-451-3943

- or -

Gary Zabow, Ph.D.
Electromagnetics Division, NIST-Boulder zabow@boulder.nist.gov
303-497-4657