

Validation of Forensic DNA Technologies:


Impacts of Thresholds on the Interpretation of Low-Template Mixtures

Catherine M. Grgicak

August 6, 2014
Boston, MA

BOSTON UNIVERSITY


NIST




Principles Behind Thresholds	
Thresholds (example values)	Principles Behind (if properly set based on lab- & kit-specific empirical data)
Analytical Threshold (e.g. 50 RFU)	Below this value, observed peaks cannot be reliably distinguished from noise
Limit of Linearity (e.g. 5000 RFU)	Above this value, the CCD can become saturated and peaks may not accurately reflect relative signal quantities (e.g., flat-topped peaks) and lead to pull-up/bleed-through between dye color channels
Stochastic Threshold (e.g. 250 RFU)	Above this peak height value, it is reasonable to assume that allelic dropout of a sister allele of a heterozygote has not occurred at that locus; single alleles above this value is single-source samples are assumed homozygous
Stutter Threshold (e.g. 15%)	Below this value, a peak in the reverse (or forward) stutter position can be designated as a stutter artifact with single-source samples or some mixtures (often higher with lower DNA amounts)
Peak Height Ratio Threshold (e.g. 60%)	Above this value, two heterozygous alleles can be grouped as a possible genotype (often lower with lower DNA amounts)
Major/Minor Ratio (e.g. 4:1)	When the ratio of contributors is closer than this value in a two-person mixture, it becomes challenging and often impossible to correctly associate genotype combinations to either the major or minor contributor


Type I Error – False labeling of noise

Analytical Threshold: Peak or Noise?



Analytical Thresholds are **applied** to ensure peaks are originating from sources other than **noise**, they are **determined** via close examination of **noise** during (and possibly after) validation





What is our goal of validating an AT?

Analytical Thresholds are **applied** to

1) ensure peaks are originating from sources other than noise and,

2) minimize unnecessary allele/locus drop-out

(Stochastic thresholds, PHR thresholds and Pr(D)'s can handle the rest)!

They are **determined** via close examination of noise during (and possibly after) validation

Analytical Threshold: How is it determined?

Use data from negatives (i.e. samples with no DNA)

- Method 1.
 - Kaiser (IUPAC 1976)
 - Long & Winefordner 1983 and Krane 2007
- Method 2.
 - Currie (IUPAC 1995)
 - Long & Winefordner 1983
- Method 3.
 - Example in SWGDAM Guidelines
- Method 4.
 - Percentile Rank

Use data from DNA dilution series

Method 5.

- Miller & Miller. *Statistics for Analytical Chemistry* (Ellis Horwood & Prentice Hall)
- IUPAC 1997 ElectroAnalytical Committee

Method 6.

- 1997 IUPAC ElectroAnalytical Committee Recommendations

Method 7.

- Performance Measurements (i.e., ROC Analysis)
- 2012 Rakay et. al. FSI: Genetics

Analytical Threshold: Color Dependent?

-Negative sample run with an internal size standard (not shown) using manufacturer's recommended protocol


Negative = extraction or amplification negative


Baseline is never below 0 RFU Processed data!


Green and Blue channels seem "quieter" than yellow and red

NIST DNA Analyst Webinar Series:
Validation Concepts and Resources- Part 1

3







Summary of Results			
Method	Validation Sample Type	Analytical Threshold for Green 5s injection example	Assumptions
1	Negatives	7	Noise is Gaussian
2	Negatives	4	Noise is Gaussian
3	Negatives	18	Noise is Gaussian
4	Negatives	6	None
5	DNA Series	31	Signal is linear wrt template and Gaussian
6	DNA Series	39	Signal is linear wrt template and Gaussian
7	DNA Series (Low-Template)	13	None
7	DNA Series (All-Templates)	50	

Analytical Threshold: Conclusions	
Baseline Noise is; <ul style="list-style-type: none">-Color dependent (low-templates)-Target dependent (high-templates)-Described by a log-normal	
To determine an AT; <ul style="list-style-type: none">-Amplify numerous samples at all targets typically seen in casework (do not only use negatives or blanks)-Analyze at 1 RFU and determine the AT-Method 7 has become the method of choice for us<ul style="list-style-type: none">-No distribution assumption-Directly measures performance of different ATs and gives the proportion of false positives AND negatives	

What is the purpose of an ST?	
<p>Stochastic Thresholds are applied to</p> <ol style="list-style-type: none">1) help us determine whether we can infer a homozygous genotype when we observe 1 peak2) minimize unnecessary false heterozygotes (i.e. too many 2p's). <p>They are determined via close examination of allele drop-out during (and possibly after) validation</p>	

Stochastic Threshold – Why do we need this?

20 µL of extract

10 µL goes into PCR

$C_n = C_0(1+E)^n$

Efficiency is not always 1. It is 1 +/- error

RFU

Cycle Number

CE Preparation & Injection Variation

Stochastic Threshold - Max Height

Allele A 179 Allele B 179

Perfect World G=A,B

Allele A 29 Allele B 179

Real World G=B,B? or B,O?

- Use the same samples used to determine AT. Amplify numerous samples at all relevant targets

Peak Height Lower RFU Allele

Peak Height Higher RFU Allele

Stochastic Threshold - Method 2 (Pr(D))

Dropout probability as a function of present allele height

Dropout probability

height of present allele

Fig. A.8. Probability dropout as a function of present allele height P(D|H).

120 150 160

Gill et al. FSI Genetics, 2009, 3, 104-111.

Methods 1 and 2 minimizes the chance of wrongly deciding a heterozygous locus is homozygous

Stochastic Threshold - Method 3

Minimizing the error of incorrectly concluding a heterozygous locus is homozygous

Minimizing the error of incorrectly concluding a homozygous locus is potentially heterozygous

At various STs, for all heterozygous loci, determine proportion of heterozygous loci falsely labeled as homozygous
For all homozygous loci, determine the proportion of homozygous loci falsely considered possible heterozygotes. Plot the proportions against each other.

AT = 30

ST=150

ST=100

151 19 151

Stochastic Threshold: Conclusions

To determine an ST;


- Amplify numerous samples at all targets typically seen in casework (can be same set you used to determine the AT)
- Analyze at your AT (not AT=1) and determine the ST
- Method 3 has become the method of choice at BU
 - Unlike ATs, STs do not seem to substantially change when determined via 3 different methods
 - No distribution assumption
 - Directly measures performance of STs and gives the proportion of false positives AND negatives


Peak Height Ratio Thresholds


Evaluate PHRs at various DNA **template levels** (e.g., dilution series of DNA). Can use the same data set used to determine AT and ST.


STs and PHRs are related.


Different PHR expectations at **different peak height ranges** may be established.












Stutter Threshold from manufacturer's validation data – locus dependent

	Locus	Stutter Threshold
<ul style="list-style-type: none">- Dependent on the length of the allele- Variance of stutter ratios may increase with target- Use the same data set used to determine PHR, ST and AT to establish stutter thresholds	CSF1PO	9.2%
	D2S1338	11.1%
	D3S1358	10.7%
	D5S818	6.8%
	D7S820	8.2%
	D8S1179	8.2%
	D13S317	8.0%
	D16S539	10.4%
	D18S51	17.0%
	D19S433	13.3%
	D21S11	9.4%
	FGA	14.7%
	TH01	5.1%
	TPOX	4.8%
	vWA	12.6%



D21S11

AT	30RFU
ST	150RFU
Stutter Filter	9.4%
PHR	0.2 (<500RFU) 0.5 (>500RFU)
Major:Minor	4:1

28 100 30 876 31 883

1. No. of contributors? 2.
2. Is the 29 an allele? Yes.
2. Possible non-observed allele if 2 contributors? No.

List out all alleles:
28, 29, 30, 31 and consider all possible pairs – assuming 2 contributors

Person 1	Person 2
28,29	30,31
28,30	29,31
28,31	29,30
29,30	28,31
29,31	28,30
30,31	28,29

D21S11

AT	30RFU
ST	150RFU
Stutter Filter	9.4%
PHR	0.2 (<500RFU) 0.5 (>500RFU)
Major:Minor	4:1

28 100 30 876 31 883

$PH_{28,29} = 100(0.2) \ln \frac{100}{0.2}$
 $PH_{28,29} = 20 \text{ to } 500 \text{ RFU}$
Since $PH_{28,29} = 120$ is btw 20-500, then 28,29 is possible combination

$PH_{28,30} = 100(0.2) \ln \frac{100}{0.2}$
 $PH_{28,30} = 20 \text{ to } 500 \text{ RFU}$
Since $PH_{28,30} = 876$ is NOT btw 20-500, then 28,30 is NOT a possible combination

$\frac{100 + 120}{100 + 120 + 876 + 883} = 0.11$

Major	Minor
30,31	28,29

Person 1	Person 2
28,29	30,31
28,30	29,31
28,31	29,30
29,30	28,31
29,31	28,30
30,31	28,29

D21S11

215

(24-600)

29,120,13,7%

(20-500)

28

100

(438-1752)

30

276

(441-1766)

31

883

AT

100RFU

ST

500RFU

Shutter Filter

9.4%

PHR

0.2 (<500RFU)

0.5 (>500RFU)

Major:Minor

4:1

$PH_{28,30} = 100(0.2) \cdot \frac{100}{0.2}$

$PH_{28,30} = 20 \text{ to } 500 \text{ RFU}$

Since $PH_{28,30}=120$ is btw 20-500, then 28,29 is possible combination

Since $PH_{28,30}=876$ is NOT btw 20-500 Then 28,30 is NOT a possible combination

Possible genotype combinations (28,29,30,31)

Person 1	Person 2
28,29	30,31
28,30	29,31
28,31	29,30
29,30	28,31
29,31	28,30
30,31	28,29

$\frac{100 + 120}{100 + 120 + 876 + 883} = 0.11$

D21S11 Inferred Genotypes

Major	30,31
Minor	28,29

D16S539

9

8,87,6,9%

1254

11

186

13

735

12,75,10,2%

1. No. of Contributors? 2.

2. Are 8 and 12 alleles? Maybe.

2. Possible non-observed alleles if 2 contributors? No.

List out all possible alleles: 8, 9, 11, 12, 13 and consider all possible pairs -- assuming 2 contributors

AT

30RFU

ST

150RFU

Shutter Filter

10.4%

PHR

0.2 (<500RFU)

0.5 (>500RFU)

Major:Minor

4:1

Possible genotype combinations if 8 contained allele (8,9,11,13)

Person 1	Person 2
8,9	11,13
8,11	9,13
8,13	9,11
9,11	8,13
9,13	8,11
11,13	8,9

Possible genotype combinations if 12 contained allele (9,11,12,13)

Person 1	Person 2
9,11	12,13
9,12	11,13
9,13	11,12
11,12	9,13
11,13	9,12
12,13	9,11

Possible genotype combinations if 8 and 12 were stutter (9,11,13)

Person 1	Person 2
9,9	11,13
9,11	9,13 or 11,13 or 13,13
9,13	9,11 or 11,11 or 11,13
11,11	9,13
11,13	9,9 or 9,11 or 9,13
13,13	9,11

D16S539

9

8,87,6,9%

1254

11

186

13

735

12,75,10,2%

1. No. of Contributors? 2.

2. Are 8 and 12 alleles? Maybe.

2. Possible non-observed alleles if 2 contributors? No.

List out all possible alleles: 8, 9, 11, 12, 13 and consider all possible pairs -- assuming 2 contributors

AT

30RFU

ST

150RFU

Shutter Filter

10.4%

PHR

0.2 (<500RFU)

0.5 (>500RFU)

Major:Minor

4:1

Possible genotype combinations if 8 contained allele (8,9,11,13)

Person 1	Person 2
8,9	11,13
8,11	9,13
8,13	9,11
9,11	8,13
9,13	8,11
11,13	8,9

Possible genotype combinations if 12 contained allele (9,11,12,13)

Person 1	Person 2
9,11	12,13
9,12	11,13
9,13	11,12
11,12	9,13
11,13	9,12
12,13	9,11

Possible genotype combinations if 8 and 12 were stutter (9,11,13)

Person 1	Person 2
9,9	11,13
9,11	9,13 or 11,13 or 13,13
9,13	9,11 or 11,11 or 11,13
11,11	9,13
11,13	9,9 or 9,11 or 9,13
13,13	9,11

D16S11 Inferred Genotypes


Major	9,13
Minor	8,11 or 11,12 or 9,11 or 11,11 or 11,13

NIST DNA Analyst Webinar Series:

Validation Concepts and Resources- Part 1

11

D16S539



1. No. of Contributors? 2.
2. Possible non-observed alleles if 2 contributors?
Yes.
List out all possible alleles:
9, 11, 13, O and consider all possible pairs
– assuming 2 contributors

D16S11 Inferred Genotypes	
Major	9,13
Minor	9,11 or 11,11 or 11,13 or 11,O

AT	100RFU	
ST	500RFU	
Stutter Filter	10.4%	
PHR	0.2 (<500RFU)	0.5 (>500RFU)
Major/Minor	4:1	

Possible genotype combinations if 8 and 12 were stutter (9,11,13)

Person 1	Person 2
9,9	11,13
9,11	9,13 or 11,13 or 13,13
9,13	9,11 or 11,11 or 11,13
11,11	9,13
11,13	9,9 or 9,11 or 9,13
13,13	9,11

Possible genotype combinations if DO (9,11,13,O)

Person 1	Person 2
9,11	13,O
9,13	11,O
9,O	11,13
11,13	9,O
11,O	9,13
13,O	9,11

Profile 1. – Minor Genotype Possibilities

Given 30 AT Threshold and interpretation standard operating procedures, the genotypes of the minor are.....

Description	D8S1179	D21S11	D7S820	CSF1PO	D3S1358	TH01	D13S317	D16S539	D2S1338
Inferred Genotypes of Minor	13,14	28,29	9,11 or 10,11 or 11,O	9,14 or 10,14 or 11,14 or 12,14 or 14,O	14,14 or 14,15 or 14,16	9,9 or 9,9.3 or 9.3,O or O,O	9,12 or 10,12 or 12,13	8,11 or 9,11 or 11,11 or 11,12 or 11,13	15,18 or 16,18 or 18,18 or 18,19 or 18,O
Known	13,14	28,29	10,11	10,14	14,16	9,9.3	12,13	11,13	18,19

Description	VWA	TPOX	D18S51	AMEL	D5S818	FGA
Inferred Genotypes of Minor	14,17 or 15,17 or 17,17	8,8 or 8,11 or 11,11 or 8,O or 11,O or O,O	14,15 or 15,15 or 15,19 or 15,20	X,Y	11,11 or 11,12 or 11,13	20,22 or 21,22 or 22,22 or 22,23 or 22,O
Known	14,15	8,11	15,20	X,Y	11,13	22,23

O= any other allele (not observed)

Conclusions

- Thresholds should be set to minimize false positives AND false negatives.
- These thresholds can be applied in a systematic way in order to deduce genotypes of 1- and 2- person mixtures
 - In particular, peak height ratio thresholds can be used to examine possible genotype combinations
- The same validation data-set of single-source samples of known genotype can be used to establish all thresholds. If possible, use actual samples (not dilution series from concentrated stock).
- Thresholds obtained for a given method must be applied only to evidence obtained using the same method (i.e. kit, injection time, etc).

NIST DNA Analyst Webinar Series:
Validation Concepts and Resources- Part 1

12