





Molar Mass and Size Measurement of Proteins via Tandem Flow Field Fractionation-Light Scattering

Wyatt N. Vreeland, Ph.D.
Chemical Engineer
Material Measurement Laboratory
National Institute of Standards & Technology
12 Feb 2015



PART I

Field Flow Fractionation



The Basics

A Brief History



- J. Calvin Giddings (1930-1996), Prof., Univ. of Utah
- Described FFF in 1976, Science (1976): 1244-1245
- Avid outdoors man



The Ground Rules

Particles must be smaller than the length perpendicular to the elution axis

$$R_{\text{particle}} \ll \text{Thickness}_{\text{channel}}$$



Channel flow must be purely laminar ($Re \ll 1$)

What FFF Can Do?

- Separates particles on the ratio of two (pseudo) orthogonal physicochemical properties
- Can separate “delicate” or “soft” colloidal particles
- Different from chromatography because there is NO STATIONARY PHASE

Typical analytical
LC
column packing
surface area =
~900
m²




Typical analytical
FFF channel
surface area =
~0.05 m²

font sized normalized to a
log₁₀ scale of relative
surface area

PART II

Absolute Size Measurement



MultiAngle Laser Light
Scattering

UV/Vis Spectrometry
Refractive index

$$R(\theta) = K^* Mc P(\theta)$$

$$K^* = \frac{4\pi^2 n_0^2}{\lambda_0^4 N_A} \left(\frac{dn}{dc} \right)^2 \quad \frac{1}{P(\theta)} = 1 + \frac{16\pi^2}{3\lambda^2} A_2 \sin^2 \frac{\theta}{2}$$

Absolute Size Determination

R_g & D_T or R_h

Radius of Gyration

Measure the *angular dependence* of scattered light

$$R(\theta) = \frac{4\pi^2 n_0^2}{N_A \lambda_0^4} \left(\frac{dn}{dc} \right)^2 M_w P(\theta) [1 - 2A_2 M_w P(\theta)]$$

Translational Diffusion Coefficient
-Or-
Hydrodynamic Radius

Measure the *time dependent fluctuation* of scattered light

$$g(\tau) = \frac{\langle I(t) \cdot I(t + \tau) \rangle}{\langle I(t) \rangle^2} = e^{-q^2 D_T \tau} \quad q = \frac{4\pi n_0}{\lambda} \sin(\theta/2)$$



A Case Study: GroEL

Tetradecamer assembly

57 kDa monomer

798 kDa complex

Nickname: “Shake-N-Bake”



dRI (concentration)



low angle LS (noisy)

medium angle LS (high S:N ratio)

high angle LS (high S:N ratio)

A Case Study: GroEL

Tetradecamer assembly
57 kDa monomer
798 kDa complex
Nickname: “Shake-N-Bake”





$R_n = 22.4 \text{ nm}$
 $R_w = 22.3 \text{ nm}$
 $R_z = 22.2 \text{ nm}$
(large uncertainty)

$M_n = 775.1 \text{ kDa}$
 $M_w = 789.6 \text{ kDa}$
 $M_z = 804.7 \text{ kDa}$
(small uncertainty)

A Case Study: GroEL

Tetradecamer assembly
57 kDa monomer
798 kDa complex
Nickname: “Shake-N-Bake”



$$D_{t \text{ (avg)}} = 4.33e^{-8} \text{ cm}^2/\text{sec}$$
$$R_h = 8.2 \text{ nm}$$

Future Directions


Inter-comparison of AIST and NIST techniques to these nascent measurement challenges

Exploration of large experimental space:

- Storage conditions
- Lyophilization protocols
- Reconstitution techniques

Move biological measurements to a more quantitative dogma

Ideas



Thoughts



Suggestions

