

The background image shows a wide-angle sunset over a calm sea. The sky is filled with warm, orange, and yellow hues, with scattered clouds reflecting the light. In the foreground, a sandy beach is visible. Two silhouetted figures are seated on the sand, facing the ocean. One figure is closer to the left, and the other is further right, both appearing to be looking at something in their hands or on the ground.

Rapid Analysis of Synthetic Polymers and Polymer Additives Using Multi-Sample Solvent-free MALDI Target Preparation and ASAP Mass Spectrometric Approaches

Charles McEwen¹ and Sarah Trimpin²

- 1) DuPont Corporate Center for Analytical Sciences, Wilmington, DE**
- 2) University of Indiana, Bloomington, IN**

Polymers

Are complex mixtures for which physical properties such as:

- O tensile strength**
- O elasticity**
- O thermal and oxidative stability**

as well as molecular properties such as:

- O repeat unit**
- O end group or pendant group compositions**
- O branching**
- O molecular weight distribution**
- O additives composition**

are required for complete characterization.

Mass Spectrometry Tools

- **MALDI**
- **LD**
- **ESI**
- **FD**
- **SIMS**
- **GC/MS**
- **LC/MS**
- **KIDS**
- **Pyrolysis (EI or CI)**
- **ESI-IMS/MS**

Challenges for Mass Spectrometry

- **Molar mass distribution** (e.g. wide PD polymers)
- **Ionization** (e.g. non-polar, polyanion, insoluble polymers)
- **Chemical composition** (e.g. high-mass homopolymers, and lower mass co- and ter-polymers, especially with multiple end groups)
- **Secondary structure** (random, block, branched)
- **Tertiary structure** (e.g. elastomers)
- **Quantitation** (e.g. ends, additives, blends)
- **Fast, cheap and accurate**

Electrospray

Peter Chen
Sarah Trimpin




ESI MS for Studying Polymerization Catalysis

- Ethylene oligomers attached to homogeneous Ziegler-Natta Pd catalysts.
- Distributions were observed for the catalytic centers before and after chain transfer
- Obtained absolute rates of initiation, propagation, and chain transfer
- Hinderling, C.; Chen, P. *Int J Mass Spectrom*, 2000, 195/196, 377-383.

Identification of Inactive Catalyst Species and Active Polymer/Catalyst Complexes

ESI of Versipol Catalysts with Growing Polymer Chain



Accelerated Articles

Resolving Oligomers from Fully Grown Polymers with IMS–MS

Sarah Trimpin, Manolo Plasencia, Dragan Isailovic, and David E. Clemmer*

Department of Chemistry, Indiana University, Bloomington, Indiana 47405

Ion mobility and mass spectrometry techniques, combined with electrospray ionization, have been used to examine distributions of poly(ethylene glycols) (PEG) with average molecular masses of 6550 and 17 900 Da. The analysis provides information about the polymer size distributions as well as smaller oligomers existing over a wide range of charge states and sizes (i.e., $[HO(CH_2-CH_2O)_xH + nCs]^{n+}$, where x ranges from 21 to 151 and $n = 2$ to 11 for the 6550 Da sample; and, x ranges from 21 to 362 and $n = 2$ to 23 for the 17 900 Da sample). The present data show that oligomer distributions also fall into families, corresponding to much narrower size distributions for individual charge states; this dramatically simplifies data analysis. For example, we show evidence for baseline resolution of the +10 charge state of polymers. Unlike the charge-state trends reported previously

and are, therefore, of significant interest.^{1–3} Such compounds are described by indication of the monomer unit composition and a value that represents the average size or molecular weight [e.g., poly(ethylene glycol) (PEG) 6550, which corresponds to a polymeric ethylene oxide chain containing ~150 monomer units]. Historically, methods aimed at characterizing polymers (e.g., light scattering or size exclusion chromatography)⁴ only provide information about the bulk material. However, it is generally understood that even simple polymers created from well-characterized reactions may exist as a mixture of sizes commonly defined by the polydispersity index (PDI).⁴ If more complex chemistries are employed, significant chemical heterogeneity may also be present, such as variations in end-groups or within the polymer chain (e.g., copolymers).




Figure 1. Schematic diagram of the ESI-IMS-MS instrument. Ions accumulated in the source funnel (F1) are gated through ion gate G1, where they undergo mobility separation in consecutive drift regions (D1, D2). Ions can be selected on the basis of mobility at G2 and used to determine cross sections. Samples are infused into the source region by means of a Nanomate autosampler mounted at the ESI source. For more detail on ion selection and cross section determination, refer to the text.






Figure 2. Two-dimensional $tb(m/z)$ distribution of PEG 6550 doped with cesium acetate. The $tb(m/z)$ distribution is shown on the left (inset A) with the integrated mass spectrum (trace on left of inset A) and integrated drift time distribution (top trace on inset A). The groups in the distribution corresponding to oligomers and polymers are labeled. Inset B (right panel) shows an expanded view of the $tb(m/z)$ region of the most intense polymer features with representative integrated drift time distribution and mass spectrum. The charge-state resolved components are labeled by charge-state families ranging from low-intensity +2 oligomers to the more intense +10 and +11 polymers (white arrows).









Figure 4. Two-dimensional $t_0(m/z)$ distribution of PEG 17 900 doped with cesium chloride. The $t_0(m/z)$ distribution is shown on the left (inset A) with corresponding total drift time distributions and integrated mass spectrum (top and leftmost trace, respectively). Prominent features in the distribution include cesium chloride clusters and regions dominated by oligomers and polymers (C, see Figure 5). Inset B shows an expanded view of the $t_0(m/z)$ region where the representative polymer (high charge-state) family is separated from smaller species (oligomers).



MALDI MS of a Random Copolymer



clemmer@indiana.edu

Matrix-Assisted Laser Desorption/Ionization (MALDI)

sample
holder



- analyte
- matrix

➡ desorption and
ionization of the
macromolecules

Properties and Significance of MALDI

- Mild desorption/ionization: often fragment-free
 - Nearly unlimited mass range
 - Broad applicability to macromolecule characterization:
 - Low sample amounts required
→ important for biological applications
 - Singly charged molecule ions
→ important for synthetic polymers and mixtures
- ⇒ MALDI and ESI complementary**

Solvent-based Sample Preparation

homogenization

transfer



dissolved
matrix

dissolved
analyte

dissolved
analyte/matrix

homogenous solution!?

homogenous crystal formation!



Solvent-free Sample Preparation

homogenization
dry matrix


dry analyte

mechanical mixing

The diagram illustrates the process of mechanical mixing. It shows two separate test tubes: one containing a yellow "dry matrix" and another containing a blue "dry analyte". An arrow points from each tube towards a central point where they converge. A pink curved arrow labeled "mechanical mixing" indicates the motion of the tubes as they come together. Below this convergence point, a green test tube contains a uniform green "dry analyte/matrix" mixture, representing the resulting homogenous powder.

homogenous
powder!?

transfer



no changes in
homogeneity!!

Analyte/matrix homogenization



Mortar and pestle




Ball Mill (mini ball mill)



Vortex devices

Advantages Beyond Solubility of Synthetic Macromolecules

- Low cost
- Green chemistry
- No side products
- No side reactions
- Oxidative stability



(stabilization)
s)
(materials)
e.g.
oxidative


Advantages Beyond Solubility of Synthetic Macromolecules

- Less laser power (fragmentation intense analytes)
- No segregation (fewer optimization trials)
- No alteration of sample during sample prep

Trimpin, S., Keune, S., Räder, H.J., Müllen, K. *J. Amer. Mass Spectrom.* 17: 661-671, 2006.


Solvent-free single sample preparation

Homogenization + transfer




Solvent-free multisample preparation

homogenization/transfer simultaneously




- Sample preparation reduced to a single step.
- Time savings
- Multiplexed sample preparation gives prospect for high throughput

Bactiplate (A) and Mirror Image MALDI Plate (B) Employing 5 min. Vortexing




Multisample Solvent free MALDI Preparation

A)




PEG 1400

B)




PMMA 1830

C)



PMMA 1830, 2400, 2990,
3800, 5270, 6950, 10300
Equi-weight

Multisample Solvent-free MALDI Preparation: IA and IB comparing reproducibility II A and II B comparing solvent-free vs. solvent-based



**PS 7200
Solvent-free**


**PS 7200
Solvent-based**

**PMMA 10300
Solvent free**


**PMMA 10300
Solvent based**

PMMA 21000

5 min
Solvent
free




Solvent
based






Multisample(25) Solvent-free MALDI Analysis

2 min




5 min



Homogeneous Surface Important for Quantitation


Solvent-free multisample preparation

homogenization/transfer simultaneously




- Sample preparation reduced to a single step.
- Time savings
- Multiplexed sample preparation gives prospect for high throughput
- Miniaturization aiding sensitivity increase

Multisample Solvent-free MALDI Analysis



custom-made sample
holder



► MALDI plate

Precise transfer

LC 100-Sample Solvent-free MALDI Analysis

Model Polymer Mixture of 5 PEG samples

Dr. Jana Falkenhagen, Dr. Steffen M. Weidner; BAM, Berlin, Germany.

LC fractionation of PEG samples



Dr. Sarah Trimpin, OSU, Corvallis, OR.


100-sample solvent-free MALDI sample preparation



Dr. Charles N. McEwen, DuPont, Wilmington, DE.

Mass Analysis

Liquid Adsorption Chromatography at “Critical’ Conditions of Adsorption (LACCC)-Fractionation of PEG samples¹



System:

Column: 2 X YMC RP 18 (150, 300 Å)
Solvent: Methanol / Water; 83/17, v/v
Detection: ELSD, SEDEX 45 (Sedere, ERC)
Temperature: 45°C
Flow: 0.5 ml/min
Concentration: 5 samples,
total amount: **5.41 mg/mL**



Injection amount: **10 μl and 1 μl**
1 fraction (25 μL) per 3 seconds

► **2 x 43 fractions collected in 200 μL vials with 3 beads present,**
dried to completeness, sent out for solvent-free multi-sample analysis

¹ Falkenhagen, J., Weidner, S.M. *Rapid Commun. Mass Spectrom.* 19: 3724-3730, 2005.

Solvent-free Multisample Loading


Automated loading




- Simplifies and automates sample loading
- Time savings
- Non ESI-sprayable LC conditions possible
- Non traditional solvent systems may aid improved separation
- Ionization of molecules inaccessible by traditional MALDI

Multiple MALDI plates prepared from the same assembly using the 100-sample homogenization/transfer method (vortex-based)

MALDI plate 1: 10 minutes




MALDI plate 2: 1 minute





Mass spectra of the MALDI plate 2 (1 minute transfer)

10 μ L injection



Peak 1


PEG 6000



Peak 4


4.3E+4

PEG1000





Mass spectra of the MALDI plate 2 (1 minute transfer)


1 μ L injection



PEG 500



PEG 1000




Conclusion: Multisample solvent-free preparation and transfer

- Total time per sample is < 1 min.
- Eliminates need to dissolve analyte
- Best for polymers with $M_n < 10 \text{ kDa}$
- Use with multiple samples or for optimization
- Prepare >100 samples simultaneous
- Smaller mixing volume translates to lower sample requirement
- All the advantages of the solvent-free method

Polymer Identification and Analysis of Additives

Atmospheric Solids Analysis Probe
Mass Spectrometry
(Fast/inexpensive)

ASAP-MS




Direct Analysis using ASAP

- **Rapid analysis of**
 - Solids
 - Liquids
 - **Polymers**
 - Biological tissue
- **Sensitive, ease of use, broad applications**
- **No interference with standard ion source operations**
- **Enclosed ion source for safety**


ASAP for Polymer Additive and Volatile Oligomer Analysis

- Obtain mass of additives in seconds
- No sample extraction or preparation required
- No vacuum lock
- Accurate mass or MS/MS for confirmation
- Molecular distribution can be obtained for volatile polymers
- No interference with ESI/APCI operation

ASAP Ion Source




Commercial ASAP Probe






ASAP Waters Fishbowl Probe

**Thermo-Fisher
Ion-Max Source
With ASAP Probe**



Nitrogen Purge Ion Molecule Reactions







Polyethylene glycol
by ASAP

1700

Fluorinated acrylate
By ASAP-MS
Failed by MALDI and ESI




F containing polymer by ASAP



Polymer Additive Analysis


- Obtain mass of additives in seconds
- No sample extraction or preparation required
- Accurate mass or MS/MS for confirmation



321.3152
C₂₂H₄₁O₁


MSMS of m/z 338 ion of Erucamide

[MH - NH₃]⁺




RT: 0.00 - 0.43


Bromonated Flame Retardant: 1 μ L of 100 ppb solution (300 fm)




12 Br



Brominated Flame retardant




ASAP MS disp. Pipette tip 300 C




ASAP of Nonwoven Fiber (250 C)

nonwoven_fabric_060118155342 #63 RT: 1.89 AV: 1 NL: 1.00E7
T: FTMS + p APCI corona Full ms [100.00-2000.00]




1. Erucamide
2. Fragment of 4
3. Irganox 1076
4. Irgafos 168
5. Irgafos 168 oxidized
6. Irganox 3114

Lab coat using ASAP MS 425C




ASAP 500 C of PET Pellet



Atmospheric Pressure Pyrolysis Ionization

ASAP on steroids

ASAP Ion Source



API pyrolysis of PET pellet at 450 C (first scan...initial temp)

577.1333
C₃₀H₂₅O₁₂
trimer

Phthalate
149.0229
C₈H₅O₃

monomer
193.0493
C₁₀H₉O₄

79.0535
C₆H₇
105.0329
C₇H₅O₁

dimer
385.0915
C₂₀H₁₇O₈
341.0655
C₁₈H₁₃O₇

219.0651
C₁₂H₁₁O₄
297.0754
C₁₇H₁₃O₅

429.1175
C₂₂H₂₁O₉

533.1072
C₂₈H₂₁O₁₁


489.1172
C₂₇H₂₁O₉

621.1591
C₃₂H₂₉O₁₃
682.1624
C₄₄H₂₆O₈


tetramer
769.1758
C₄₀H₃₃O₁₆

m/z


ASAP 500 C of PET Pellet




API pyrolysis of PET pellet 450C




Lab Coat by API Pyrolysis 700C




Carpet pyrolysis API 425 C




Carpet by ASAP 500 C




Pyrolysis API of carpet 700 C




Pipette tip API pyrolysis 700 C




API Pyrolysis pipette tip (polypropylene)




Polystyrene standard 35K by ASAP 450 C



API Pyrolysis PS Standard 425 C



Crude Oil Sample by ASAP MS



Conclusion: ASAP and Py-ASAP MS of Polymers

- **Rapid methods for analysis of additives and identifying polymers**
- **AP pyrolysis gives better representation of high mass species than traditional pyrolysis methods**
- **Useful in combination with TGA analyses**
- **Unknown: What ions will be observed below m/z 50**

ASAP for Polymer Additive and Volatile Oligomer Analysis

- Obtain mass of additives in seconds
- No sample extraction or preparation required
- No vacuum lock
- Accurate mass or MS/MS for confirmation
- Molecular distribution can be obtained for volatile polymers

Acknowledgements

Richard G. McKay
Barbara S. Larsen
Patricia M. Peacock

DuPont Corporate Center for Analytical Sciences